# Weight sorting both tablets and capsules

# **CASE STUDY**



#### **SUMMARY**

Almac's SADE SP440 is used in R&D applications for weight sorting both tablets and capsules, where batch sizes can vary from 1,000-250.000

#### **CUSTOMER**

Almac Pharma Services provides the most comprehensive range of services extending from research through pharmaceutical and clinical development to commercialization of product.

#### **CHALLENGE**

- Versatile to handle both tablets and capsules without change parts
- Flexible performance for batch sizes that varies from 1,000 to 250,000

### SOLUTION

Offering high accuracy and the ability to handle many shapes and sizes without change parts, the SADE SP440 is ideal for the 100% check-weighing of large scale R&D batches.

## **OUTCOME**

"SADE equipment is used by our Clinical Services Division. They recommended CI Precision to us and we now have a SADE SP140 and a SP440 in Pharma Services"

reliable and suitable for our applications





weight sorting